



biocrates
The future of research and health

Standardized metabolomics kits for quantitative and reproducible results


Kit workflow in the laboratory and
the MetIDQ™ software

biocrates.com

Kit workflow and the MetIDQ™ software

Content

- Company overview
- Targeted metabolomics
- Standardized assay technology
- AbsoluteIDQ® p180 kit I
 - Key benefits
 - Analytical coverage and performance
 - How it works
- AbsoluteIDQ® p180 kit II
 - Setup
 - Workflow
 - Assay preparation
- MetIDQ™ software – live demo



Company overview

Company profile
Technology

Company profile



Biotech company for targeted metabolomics kit technology in Innsbruck, Austria



> 15 years of experience, biocrates technology was applied for thousands of publications and by hundreds of clients in pharma, academia and health institutions



Fully integrated mass spectrometry-based solutions for metabolic & functional microbial phenotyping across all indication areas




Kit technology empowers quantitative analysis of thousands of endogenous metabolites. Statistical analysis and biological interpretation is available.



Our customer base includes global tier one pharma and academic customer as well as national health organizations.
A global community of more than 100 core facilities, CRO's, service customers & kit users to accelerate your science

Technology

biocrates develops, produces and distributes robust and high-performance kit tools for research purposes. This is THE standard for future diagnostic applications.



Targeted metabolomics


Metabolomics research

Why metabolomics matters

Metabolomic data interpretation

Metabolomics research

Metabolomics is the analysis of the transition of small molecules along defined pathways




Typical metabolites

- Sugars
- Nucleosides
- Organic acids
- Ketones
- Aldehydes
- Amines

- Amino acids
- Small peptides
- Lipids
- Sterols
- Terpenes
- Alkaloids

- Metabolomics is the study/analysis of metabolites in biological systems
- Metabolites are low molecular weight organic molecules
- Metabolome is the total metabolites pool




Why target metabolomics matters

Our approach

- Metabolomics is the most relevant tool to answer the multidimensional questions of system biology.
- From just a few microliters of sample material, the phenotype of an individual can be uncovered.
- Quantitative, standardized, and quality-controlled targeted metabolomics methods guarantee robust and reproducible results that are the cornerstones of all clinical research.
- biocrates provides mass spectrometry-based targeted metabolite analysis, with assays covering over 1,000 metabolites from more than 30 metabolite classes.


Why metabolomics matters – The omics map

Metabolomics is functional readout of integrated genomic information under impact from exposome



Why metabolomics matters – The source

Liquid biopsies provide access to pathophysiological status of individuals




Systemic metabolism reflects energy homeostasis, immune response to disease and also functional impact of the

Metabolomics data interpretation

Disease know-how drives “metabolic signature”

Combined statistics and biochemical expertise required for state-of-the-art analysis



Standardized assay technology

Technological benefits


We standardized metabolomics

Scientific impact

Metabolomics benchmark

Technological benefits

Connect the dots to bring research to a new level – Measure locally and connect science globally



Total ownership of the results within a global network – Realizing the FAIR principles

Standardized targeted metabolomics

Standardized metabolomic data connect science globally



Kit components

- Reagents & consumables
- Methods & protocols
- Workflow manager MetIDQ

Laboratory environment


- Inter-laboratory



- Inter-MS instrument



- All different technologies delivering absolute concentrations



Database MetIDQ

Harmonized & standardized data format in one data base

long time reliable for re-analysis



Combined analysis with other omics data (e.g. GWAS)

Scientific impact

> 1,000 publications, > 27,000 citations and h-Index 79 support strong scientific impact of our technology



| | | | |
|------|------------------------------|----------------------------|--------------|
| 2018 | European Heart Journal | Tomas et al. | Cardiology |
| 2018 | Cell Metabolism | Della Torre et al. | Hepatology |
| 2018 | Cell Metabolism | Raud et al. | Immunology |
| 2018 | Gut | Mayerle et al. | Oncology |
| 2018 | Alzheimer's & Dementia | Mahmoudian Dehkordi et al. | Neuroscience |
| 2018 | Alzheimer's & Dementia | Tynkkynen et al. | Neuroscience |
| 2018 | J Cachexia Sarcopenia Muscle | Pin et al. | Oncology |
| 2018 | Nature Communications | Mao et al. | Aging |

Impact factor
> 12.000



| | | | |
|------|----------------------------------|---------------|-----------------|
| 2018 | PLoS Medicine | Varma et al. | Neuroscience |
| 2018 | JNCI | Lu et al. | Oncology |
| 2018 | PNAS | Willis et al. | Systems Biology |
| 2018 | PNAS | Skene et al. | Lifestyle |
| 2017 | Journal of Experimental Medicine | Bakiri et al. | Oncology |
| 2017 | European Respiratory Journal | Naz et al. | Pulmonology |

Impact factor
> 9.500



| | | | |
|------|---------------------------------|----------------------|-----------------|
| 2018 | EMBO Molecular Medicine | Jia et al. | Pulmonology |
| 2018 | PLoS Biology | Diessler et al. | Somnology |
| 2018 | Cancer Research | Mycielska et al. | Oncology |
| 2018 | BMC Medicine | DiBattista et al. | Epidemiology |
| 2018 | BMC Medicine | Lau et al. | Epidemiology |
| 2018 | Am J of Human Genetics | Romano et al. | Systems Biology |
| 2018 | EMBO Reports | Audano et al. | Systems Biology |
| 2018 | JASN | Li et al. | Nephrology |
| 2018 | Theranostics | Schaarschmidt et al. | Hepatology |
| 2018 | Int Journal of Epidemiology | Iqbal et al. | Epidemiology |
| 2018 | Cell Reports | Walters et al. | Aging |
| 2018 | International Journal of Cancer | Yuan et al. | Oncology |

Impact factor
> 7.000

AbsoluteIDQ® p180 kit I

Key benefits

Analytical coverage


How it works

Analytical performance

Matrices

p180 kit – Key benefits

Standardized quantitative metabolomics – Up to 188 quantitative metabolic features




p180 kit – Key benefits



All in one box



Excellent longitudinal & inter-laboratory reproducibility



High throughput & small sample volume of 10 µL



Standardized, simple & automated workflow



Kit components:

- Reagents & Consumables
 - Kit plate
 - Standardized samples: internal standards, calibration standards, and QCs
- Instrument methods
 - Ready to use LC-MS methods
 - **Method optimization not required**
- Documentation and support
 - Detailed user manuals
 - Video tutorials
 - Application notes for multiple sample matrices
 - **biocrates support by specialists**
- Workflow manager MetIDQ™ software
 - Project, sample, and plate management
 - Automated quantification and quality assessment
 - Data normalization
 - Export of results for easy data analysis

p180 kit – Analytical coverage

188 metabolites | lipids from 7 analyte classes

LC part

42 small molecules (2 classes)

- Amino acids (21)
- Biogenic amines (21)



analytical column

FIA part

146 lipids & hexose (5 classes)

- Acylcarnitines (40)
- Phosphatidylcholines (76)
- Lysophosphatidylcholines (14)
- Sphingomyelins (15)
- Monosaccharide (1)



Quant 500 kit – Analytical coverage

630 metabolites/lipids from 26 analyte classes

Nutrition-microbiome-host interaction

107 small molecules (14 classes)

- Alkaloids (1)
- Amine oxides (1)
- Amino acids (20)
- Amino acid related (30)
- Bile acids (14)
- Biogenic amines (9)
- Carbohydrates and related (1)
- Carboxylic acids (7)
- Cresols (1)
- Fatty acids (12)
- Hormones and related (4)
- Indoles and derivatives (4)
- Nucleobases and related (2)
- Vitamins and cofactors (1)

produced or modified by gut microbiota

523 lipids (12 classes)

- Acylcarnitines (40)
- Phosphatidylcholines (76)
- Lysophosphatidylcholines (14)
- Sphingomyelins (15)
- Cholesteryl esters (22)
- Ceramides (28)
- Dihydroceramides (8)
- Hexosylceramides (19)
- Dihexosylceramides (9)
- Trihexosylceramides (6)
- Diglycerides (44)
- Triglycerides (242)


+ 234 Sums & Ratios

MetaboINDICATOR



p180 kit – Analytical coverage – Pathways


Broad coverage of metabolic pathways



p180 kit


p180 kit – How it works, the workflow

80 samples in 22 hours



Quant 500 kit – How it works, the workflow

80 samples in 42 hours




p180 kit – Analytical performance

Ring trial across 6 international laboratories

Inter-laboratory reproducibility of a targeted metabolomics platform for analysis of human serum and plasma

Alexandros P. Siskos,^{1§} Pooja Jain,^{1§} Werner Römisch-Margl,² Mark Bennett,³ David Achaintre,⁴ Yasmin Asad,⁵ Luke Marney,⁶ Larissa Richardson,⁶ Albert Koulman,⁶ Julian L. Griffin,⁶ Florence Raynaud,⁵ Augustin Scalbert,⁴ Jerzy Adamski,^{7,8,9} Cornelia Prehn,⁷ Hector C. Keun^{1*}

Inter-Laboratory Precision (CV in %)



Number of laboratories: 6

- Imperial College, UK
- Cambridge University, UK
- Institute of Cancer Research, UK
- Intern. Agency for Research on Cancer, France
- Helmholtz Zentrum München, Germany
- Biocrates, Innsbruck, Austria


Inter-laboratory performance:

CV 7.6 %

Anal. Chem. 2017, 89, 656–665

p180 kit – Analytical performance

Proven analytical accuracy



NIST SRM 1950

“

“The field **MUST become more quantitative** if findings are to be translated to practical, clinical applications.”

David Wishart (2015)

”

“

“If we could get over this barrier [to not quantitate our data] then it will make **exchanging data between labs and studies** so much easier, and we could really move forward as a community.”

Jules Griffin (MetaboNEWS, Sept 2017)


”

p180 kit – Matrices and species

Connecting biological systems

Matrices

- Plasma, serum (only 10 µL)
- Dried blood spots
- Cell culture medium
- Tissue
- Tumor tissue:
- Lung lavage (BALF)
- Skin samples
- Blister liquid / skins
- Feces
- Follicular fluid
- Milk
- Urine
- CSF
- Saliva
- Cells
- Cell culture supplements




Species

- Human
- Mouse
- Rat
- Monkey
- Cow
- Sheep
- Pig
- Dog
- Chicken
- Horse
- Rabbit
- Zebrafish
- Chinese hamster
- C. elegans
- Soy
- Yeast: *Pichia pastoris*
- E. coli

p180 kit – Matrices and species

Connecting biological systems



AbsoluteIDQ® p180 kit II

Setup


Workflow – preparation, measurements, data processing

Assay preparation – hints and tricks

p180 kit – Workflow

Setup – first time kit use only:

I. LC-MS instrument method setup



Method setup:

- ✓ Instrument methods included in kit
No tuning or optimization required
- ✓ Install methods according to user manual
- ✓ Assistance by **customer support** on request

p180 kit – Workflow

Setup – first time kit use only:

II. System suitability test (SST)



testmix
LC



testmix
FIA

Method setup:

- ✓ Ready to use **testmix**
 - ✓ Instrument cleaning instructions
1. Inject testmix
 2. Adjust retention times (RTs)
 3. Check instrument performance (SST)
 4. Performance assessment together with **customer support**
 5. Run the starter kit

p180 kit – Workflow


1. Plate layout

2. SST & kit preparation



3. LC-/FIA-MS run




30 min



4 h



15.5 h



2 h



Total time: 22 h

6.
Interpretation

5. Performance assessment
& export of results

4. MS data processing




MetIDQ
Oxygen

p180 kit – Workflow

1. Plate layout

Kit plate design:

































































































- All kits use 96 well plates
- Some wells are used for kit specific samples
 - blank
 - zero
 - calibration standards
 - quality controls (QCs)











p180 kit – Workflow

1. Plate layout

Kit plate design:

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|---|---|---|---|---|---|---|---|---|--|---|---|---|
| A |  |  |  |  |  |  |  |  |  |  |  |  |
| B |  |  |  |  |  |  |  |  |  |  |  |  |
| C |  |  |  |  |  |  |  |  |  |  |  |  |
| D |  |  |  |  |  |  |  |  |  |  |  |  |
| E |  |  |  |  |  |  |  |  |  |  |  |  |
| F |  |  |  |  |  |  |  |  |  |  |  |  |
| G |  |  |  |  |  |  |  |  |  |  |  |  |
| H |  |  |  |  |  |  |  |  |  |  |  |  |

-  blank
 -  zero samples
 -  calibration standards
 -  internal standards
 -  up to 82 samples
- QCs in 3 levels:
-  low
 -  medium
 -  high


78 – 82 samples can be analyzed per kit

p180 kit – Workflow

1. Plate layout

Recommended: 5x QCs replicates

Best normalization performance = best longitudinal comparability



78 samples can be analyzed per kit

p180 kit – Workflow



1. Plate layout

5x QCs replicates required for:

- Highly robust, quantitative results
- Excellent accuracy & precision



- Outstanding inter-instrument, inter-laboratory, long-term reproducibility
- Applicable to multi-center, longitudinal epidemiological studies



Siskos et al., Anal Chem 2017; 89(1):656-65


p180 kit – Workflow

1. Plate layout



1. Register samples
2. Create plate layout
3. Export layout for kit pipetting and LC-MS runs

MetIDQ™ software



p180 kit – Workflow

2. SST & kit preparation

System suitability test (SST)



testmix
LC




testmix
FIA

- ✓ If required, follow instrument cleaning instructions
- ✓ Use **testmix**
 1. Inject testmix
 2. Check retention times (RTs)
 3. Check instrument performance (SST)
 4. Run the kit

p180 kit – Workflow

2. SST & kit preparation

Kit plate and components




holder with FDA-approved filter

filter

upper plate
sample preparation

capture plate
extracts for
measurements



Advantages:

- Internal standards on plate
- Easy and fast sample preparation
- Efficient derivatization, drying and extraction
- Precipitated proteins are filtered



- Calibration standards (7 levels)
- Quality controls (3 levels)
- Testmix (2 vials)

p180 kit – Workflow


2. SST & kit preparation

Kit preparation



4 h

- 10 μ L sample
- 10 μ L ISTDs
- 30 min N₂ drying
- PITC derivatization
- 25 min incubation
- 60 min N₂ drying
- Extraction
- 30 min shaking
- Filtration
- Centrifugation or positive pressure
- Split extract
- Transfer to empty plates and dilute separately




p180 kit – Workflow

3. LC-/FIA-MS run




MS
methods
included



LC-MS

- 6 min / sample
- 42 analytes
amino acids, biogenic amines




FIA-MS

- 3 min / sample
- 146 analytes
acylcarnitines, hexose, glycol-&
lysophospholipids, sphingomyelins

p180 kit – Workflow


4. MS data processing



p180 kit – Workflow

4. MS data processing

Quantification uses internal standards (ISTDs)
= “automated” intra-plate normalization



p180 kit – Workflow


4. MS data processing

- LC part:
- 7-point calibration + internal standards
 - Internal standards → intra-plate normalization
 - External standards (pure metabolites) → 7-point calibration
 - Quality controls → performance assessment + inter-plate normalization

| Analyte Selection | | Curve Type | | Weighting Type | |
|-----------------------|--|--|--|--|--|
| Name | | <input type="radio"/> LINEAR | | <input type="radio"/> No Weights | |
| Glu | | <input type="radio"/> LINEAR ZERO | | <input checked="" type="radio"/> 1/X | |
| Gly | | <input checked="" type="radio"/> QUADRATIC | | <input type="radio"/> 1/X ² | |
| His | | | | | |
| Ile | | | | | |
| Leu | | | | | |
| I ve | | | | | |
| Apply To All Analytes | | | | | |

| Calibration Standards | | | | | | |
|-------------------------------------|----------------|------|-------------|----------|----------------|------------------|
| Active | Sample Barcode | Well | Sample Type | Accuracy | Expected Conc. | Calibrated Conc. |
| <input checked="" type="checkbox"/> | 1020351488 | 49 | Standard L1 | 101.23 % | 25.00000 | 25.30756 |
| <input checked="" type="checkbox"/> | 1020351492 | 61 | Standard L2 | 101.09 % | 50.00000 | 50.54544 |
| <input checked="" type="checkbox"/> | 1020351503 | 73 | Standard L3 | 97.27 % | 250.00000 | 243.16843 |
| <input checked="" type="checkbox"/> | 1020351518 | 85 | Standard L4 | 99.43 % | 500.00000 | 497.13562 |
| <input checked="" type="checkbox"/> | 1020351522 | 2 | Standard L5 | 99.95 % | 1000.00000 | 999.53314 |
| <input checked="" type="checkbox"/> | 1020351537 | 14 | Standard L6 | 102.32 % | 1500.00000 | 1534.78711 |
| <input checked="" type="checkbox"/> | 1020351541 | 26 | Standard L7 | 98.72 % | 2000.00000 | 1974.32153 |


| QCs | | | | | | |
|----------------|------|-------------|----------|----------------|------------------|------|
| Sample Barcode | Well | Sample Type | Accuracy | Expected Conc. | Calibrated Conc. | Note |
| 1020356396 | 38 | QC Level 1 | 91.54 % | 147.00000 | 134.56985 | |
| 1020356407 | 50 | QC Level 2 | 88.84 % | 284.00000 | 252.31444 | |
| 1020356411 | 62 | QC Level 3 | 93.69 % | 1367.00000 | 1280.74573 | |
| 1020356407 | 87 | QC Level 2 | 83.42 % | 284.00000 | 236.92676 | |



p180 kit – Workflow

4. MS data processing

- FIA part:
- 1-point calibration + internal standards
 - Internal standards → intra-plate normalization
 - Quality controls → performance assessment + inter-plate normalization



p180 kit – Workflow

4. MS data processing

Performance:

LC part




quantitative
accurate & precise


Inter-plate data comparability:

- QC normalization **recommended**

FIA part



quantitative
accurate & precise




relative-quantitative
precise! accuracy?

Inter-plate data comparability:

- QC normalization **required**

p180 kit – Workflow

5. Performance assessment & export of results



| Thr | Trp | Tyr | Val | 1-Met-His | 3-Met-His | 5-AVA | AABA |
|------------|------------|------------|------------|----------------|----------------|----------------|----------------|
| Aminoacids | Aminoacids | Aminoacids | Aminoacids | Aminoacids ... | Aminoacids ... | Aminoacids ... | Aminoacids ... |
| 1.69 | 1.62 | 1.36 | 5.07 | 0.015 | 0.004 | 0.043 | 0.696 |
| 148 | 48.6 | 69.7 | 197 | 4.86 | 5.51 | 2.51 | 17.3 |
| 306 | 230 | 239 | 411 | 9.87 | 11.5 | 4.79 | 32.9 |
| 150 | 50.4 | 66.8 | 153 | 7.04 | 2.54 | 0.671 | 4.83 |
| 111 | 39.6 | 103 | 121 | 3.90 | 1.88 | 0.259 | 10.1 |
| 160 | 72.2 | 69.3 | 162 | 9.10 | 3.49 | 0.072 | 6.27 |
| 123 | 51.9 | 63.1 | 167 | 6.89 | 2.34 | 0.026 | 5.70 |
| 91.5 | 46.6 | 40.3 | 107 | 7.78 | 2.92 | 0.082 | 8.43 |
| 182 | 72.3 | 57.0 | 150 | 4.84 | 1.80 | 0.043 | 4.80 |
| 122 | 42.5 | 42.3 | 123 | 4.10 | 1.44 | 1.18 | 4.91 |
| 115 | 56.5 | 41.7 | 113 | 6.42 | 2.22 | 0.025 | 8.21 |
| 161 | 59.0 | 87.4 | 177 | 6.01 | 2.67 | 0.043 | 8.16 |
| 117 | 41.2 | 75.6 | 153 | 7.28 | 2.85 | 0.020 | 3.54 |

Automated technical validation:



- ✓ blank checking cleanliness of system
- ✓ sensitivity checking internal standard intensities
- ✓ accuracy checking quality control results
- ✓ monitoring inter-plate monitoring of QC results and CVs

Analytical approval of the plate
Data normalization (recommended)
Export of results table in various file formats

p180 kit – Workflow

6. Interpretation

What is a metabolism indicator?



A metabolism indicator can act as a proxy for enzyme activity
≠
Metabolism indicators allow to measure enzymatic activity with metabolomics

p180 kit – Workflow

6. Interpretation

Make data coherent – Why quantitation matters

Fischer ratio

↓ ratio

↑ severity of liver failure

BCAA (Val + Leu + Ile)

AAA (Phe + Tyr + Trp)



Fischer et al., Am J Surg 1974;127:40.


IDO activity

↑ ratio

↑ cancer invasiveness and progression

Kynurenine (Kyn)

Tryptophan (Trp)



Hornýák et al. 2018, Front Immunol 2018; 9:151.

Phospholipase A2 activity

↑ ratio

↑ inflammation and CVD

Lysophosphatidylcholines (LPCs)

Phosphatidylcholines (PCs)



Paapstel et al. Nutr Metab Cardiovasc Dis 2018; 28(1):44–52.

AAA: aromatic amino acids, BCAA: branched-chain amino acids, IDO: indoleamine 2,3-dioxygenase, CVD: cardiovascular disease

p180 kit – Workflow

Assay preparation – hints and tricks



p180 kit – Workflow

Assay preparation – hints and tricks

Kit plate - pipetting of samples:

- biocrates kit plate used
- sample volume: 10 μ L




Let's have a look at the pipetting of biocrates kit plates




p180 kit – Workflow

Assay preparation – hints and tricks




p180 kit – Workflow

Assay preparation – hints and tricks



p180 kit – Workflow

Assay preparation – hints and tricks




p180 kit – Workflow

Assay preparation – hints and tricks



p180 kit – Workflow

Assay preparation – hints and tricks



p180 kit – Workflow

Assay preparation – hints and tricks

Standard and maximum sample volumes:


- Standard sample volume: 10 μL
 - for plasma and serum
 - unless otherwise noted in SOP
- Further matrices: what is the optimal sample volume?
 - perform pilot test using a small number of samples
 - try different sample volumes, e.g. 10, 20, 30 μL
 - for assistance, customer should contact **customer support**



p180 kit – Workflow

Assay preparation – hints and tricks


After sample preparation...



p180 kit – Workflow


Assay preparation – hints and tricks

...drying of kit plate




p180 kit – summary

Lab workflow



p180 kit – summary

Measurements



Instrument methods included in kit:

- Ready to use
- No tuning or optimization required

Method setup:

- Install methods according to user manual
- Adjust retention times (RTs) using testmix
- Instrument test: perform SST
- Assistance by [customer support](#) on request


p180 kit – summary



All in one box



Excellent longitudinal & inter-laboratory reproducibility



High throughput & small sample volume of 10 µL



Standardized, simple & automated workflow



Key features:

- 10 µL sample volume
- multiple matrices per plate
- up to 80 samples per kit
- 4 h sample preparation time, 2 hours hands-on
- one day sample analysis run time
- 2 h data automatized analysis time



high throughput


p180: up to 500 samples per week

MetIDQ™ software live demo



Workflow manager for targeted metabolomics

- Before sample measurement:
 - register samples and projects,
 - generate plate layout and worklist for mass spectrometer
- After sample measurement:
 - quantification and data processing
 - export of results



MetIDQ



3 minutes break...



3 2 1

Stay tuned!

Contact

Dr. Manuel Kratzke

Specialist | Applications Support

+43 512 579823

manuel.kratzke@biocrates.com



© 2020 biocrates life sciences ag, an Austrian corporation.

All rights are reserved. The name biocrates and the logo are registered trademarks of biocrates life sciences ag.