

# Capillary Electrophoresis-Mass Spectrometry (CE-MS)

Coral Barbas  
Antonia Garcia



## Capillary Electrophoresis



- Separation based on electrophoretic mobility
- Simple instrumentation
- Primary applications in bioanalysis
  - DNA sequencing
  - DNA fragment analysis
- Multiple modes for improved selectivity of neutrals
  - MEKC
  - CEC



# Advantages and Disadvantages of CE

## Advantages

- Offers new selectivity, an alternative to HPLC
- Easy and predictable selectivity
- High separation efficiency ( $10^5$  to  $10^6$  theoretical plates)
- Small sample sizes (1-10  $\mu$ l)
- Fast separations (1 to 45 min)
- Can be automated
- Quantitation (linear)
- Different “modes” (to be discussed)


## Disadvantages

- Cannot do preparative scale separations
- Low concentrations and large volumes difficult
- “Sticky” compounds
- Species that are difficult to dissolve
- Reproducibility problems

# Applications of CZE

- Wide variety of applications
  - Small molecules
  - Macromolecules (proteins, peptides)
- Limitations
  - Must have different charges
  - Low ionic strength sample
- Advantages
  - Simple
  - Direct analysis of complex systems


# Capillary Electrophoresis (CE)



## Modes of CE


- Capillary Zone Electrophoresis (CZE)
  - Basic mode using open channels
- Micellar Electrokinetic Chromatography (MEKC)
  - Separates compounds with micelles
- Capillary Gel Electrophoresis
  - Size exclusion using sieving gels
- Capillary Electrochromatography
  - Hybrid of CE and HPLC
- Capillary Isoelectric Focusing
- Enantiomeric CE

# Capillary Zone Electrophoresis (CZE)



## Electroosmotic Flow


- Capillary flow mechanism based on applied potential and pH
- Provides bulk solution flow in capillaries with moderate to low concentration buffers are used
- Very dependent on solution ionic strength and surface chemistry



# Electrophoretic Mobility

$$\mu = \frac{q}{6\pi\eta r}$$

$\mu$  = electrophoretic mobility  
 $Q$  = charge on the particle  
 $\eta$  = solution viscosity  
 $r$  = Stokes radius of the particle




# Electrophoresis and Electroosmosis


A pictorial representation of the combined effect in a capillary, when EO is faster than EP (the common case):

$$v = (\mu_{ep} + \mu_{eo})E = (\mu_{ep} + \mu_{eo})\frac{V}{L}$$

Figure from R. N. Zare, Stanford




# Diagrams of flow in CE and HPLC



## Detection Options


Detector	Advantages	Characteristics	Limit the range of detection (M)
UV/vis absorbance	-The possibility of direct and indirect detection -Very common detector	-Universal	- $10^{-3}$ - $10^{-6}$ for the detection of aromatic compounds
LIF (Laser-induced fluorescence)	-Highly sensitive and highly selective -Used for fluorescent compounds or derivatives	-Selective	- $10^{-6}$ - $10^{-9}$
MS	-Qualitative and quantitative information -Highly sensitive and highly selective	-Universal -Selective	- $\approx 10^{-5}$ (it depends on the type of MS and metabolites)



# Optimizing CE Separations PARAMETERS


- pH
  - First parameter to control
  - Effects EOF and mobility (charge)
- Organic Solvent
  - Analyte solvation
- Interacting agent
  - Ion-pairing, solvation, etc.
- Non-aqueous Conditions
  - Solvation and charge
- Temperature
  - Solvation, chemical equilibria

## INSTRUMENTATION CE-MS



\* OPTIMIZATION REQUIRED

# ESI SOURCE CE-MS



## CE-MS: Electrical interfacing




image provided by Agilent Technologies

---

This project has been funded with support from the European Commission.

This publication reflects the views only of the authors, and the Commission cannot be held responsible for any use which may be made of the information contained therein

